

## KAINOS®+ bone graft substitute Instructions for Use

### Device description:

KAINOS+ is a microporous and macroporous biphasic calcium phosphate ceramic consisting of 60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP). KAINOS+ is available in various shapes and sizes. KAINOS+ may be used with physiological saline, patient's own serum, whole blood, or bone marrow aspirate (BMA).

KAINOS+ is provided sterile for single patient use.

### Indications for use:

KAINOS+ is intended for use as a bone void filler for bony voids or gaps of the skeletal system (e.g. extremities, spine and pelvis) that are not intrinsic to the stability of the bony structure. KAINOS+ can be used with autograft as a bone graft extender. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone.

KAINOS+ is a bone void filler without initial mechanical properties. Therefore rigid fixation techniques may often be recommended. When packed into a bony site, KAINOS+ gradually resorbs and is replaced with new bone during the healing process.

### Caution:

- KAINOS+ is intended for single use only and must not be re-used. The re-use of an implant may result in its failure, in infection and/or death.
- Federal law restricts this device to sale by or on the order of a physician.

### Contraindications:

KAINOS+ has limited initial mechanical properties. Therefore, this product is contraindicated where the device is intended as structural support in the skeletal system.

Conditions representing contraindications include also:

- Osteomyelitis
- Implantation in necrotic surgical sites
- Degenerative bone disease
- Intra-articular implantations

### Complications:

Possible adverse effects include, but are not limited to:

- Wound complications including hematoma, infection, and other complications that are possible with any surgery
- Incomplete, or lack of, osseous ingrowth into bone void, as is possible with any bone void filler

### Precautions:

- Rigid fixation techniques may be required to assure rigid stabilization of the defect in all planes.
- Maximum contact between the product and the recipient bone must be established.
- The implantation in a revision surgical site containing non-resorbable fragments of material (e.g. polyethylene ligament waste, carbon fibers) is not recommended.
- As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure. This includes (but is not limited to) individuals with longterm steroidal therapy or treatment acting on the calcium or phosphorus metabolism.
- KAINOS+ is radiopaque until resorbed. Radiopacity may mask underlying pathological conditions. Radiopacity may also make it difficult to radiographically assess the ingrowth of new bone.

### Warnings:

- The implant cannot be re-sterilized.
- Reprocessing and/ or reuse may result in infection and/ or loss of function, which ultimately may result in the death of the patient.
- KAINOS+ should only be used by surgeons experienced with orthopedic fixation devices and bone grafting.

### Material:

KAINOS+ is a biphasic calcium phosphate ceramic consisting of 60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP). The chosen materials are biocompatible and MRI-compatible.

### How supplied:

The KAINOS+ devices are provided sterile. All sterile products are labeled "STERILE." Do not use sterile products if the packaging has been damaged or previously opened. Do not re-sterilize or autoclave sterile implants.

### Preparation:

- Use care to avoid destruction of porous structure.
- KAINOS+ may be soaked in patient's serum, whole blood, or bone marrow aspirate (BMA) until graft is completely impregnated.
- Prior to this, KAINOS+ must be hydrated without excess with physiological saline to prevent osmotic damage.
- Based on experiments, the ideal proportion of fluid should be 1 volume fluid to 2 volumes of KAINOS+.
- Prepared KAINOS+ should be used immediately to preserve cell viability.

### Handling and use:

KAINOS+ is provided sterile and should be considered sterile unless the inner packaging has been opened or damaged. This product must not be re-sterilized. The sterilization process of this product using radiations was validated during its conception. We do not guarantee the sterility of the product if a new sterilisation is performed on it. Precautions must be applied in order to preserve the porous structure. This device is for single patient use only and should never be reused. The non-respect of the single use can conduct to septic issues. Additionally, the resorption process of this bone substitute started immediately after its implantation on the osseous defect so it can not be re-used.

KAINOS+ should be stored at ambient temperature.

### Further Information:

Recommended directions for the use of these systems (surgical operative techniques) are available at no charge upon request. For additional information please contact:



SIGNUS MEDICAL LLC

18888 Lake Drive East

Chanhassen, MN 55317, USA

Phone: 1-952-294-8700

Fax: 1-952-975-0465

E-mail: [info@signusmedical.com](mailto:info@signusmedical.com)

### Product Guarantee:

SIGNUS Medizintechnik GmbH guarantees that each individual spinal implant is manufactured with the greatest of care and from selected material. The entire process, from manufacture to final packaging, is under constant quality control. However, given that SIGNUS Medizintechnik GmbH has no influence over the selection or use of the implant, the diagnosis of the patient or the surgical technique used, or the handling of the implant following dispatch from our company, no guarantees can be given regarding a successful surgical result or the lack of complications.